

ProteCol™ Guard Columns

SGE's unique guard column design couples directly to the analytical column



Part no.	Particle Size (µm)	ID	Length (mm)	Pore Size (Å)	Phase Type	Description	Pack Size
250033	5	4mm	10	120	C18 GP	Guard Column	3
250003	5	2.1mm	10	120	C8 HQ	Guard Column	3
250005	5	4mm	10	120	C8 HQ	Guard Column	3
250007	5	2.1mm	10	100	C18 HQ	Guard Column	3
250009	5	4mm	10	100	C18 HQ	Guard Column	3
250013	5	2.1mm	10	120	C18 HPH	Guard Column	3
250015	5	4mm	10	120	C18 HPH	Guard Column	3
250019	3	2.1mm	10	1000	C8 HQ	Guard Column	3
250017	3	4mm	10	1000	C8 HQ	Guard Column	3
250021	3	2.1mm	10	200	C18 HQ	Guard Column	3
250023	3	4mm	10	200	C18 HQ	Guard Column	3
250025	5	2.1mm	10	300	C18 HQ	Guard Column	3
250027	5	4mm	10	300	C18 HQ	Guard Column	3
250029	3	2.1mm	10	300	C18 HQ	Guard Column	3
250031	3	4mm	10	300	C18 HQ	Guard Column	3

Phase description: HQ = High Quality, GP = General Purpose, HPH = High PH