

FMS-15 Ultra-thin magnetic stirrer

Magnetic Stirrer

Introduction

Ultra-thin magnetic stirrer allow precise and stable stirr of liquids over a wide range, especially for small volume samples, It is a necessary tool for petroleum, chemical industry, environmental protection, education and scientific research, etc.

Feature

1. Ultra-light and Ultra-thin, easy to store and save laboratory space.
2. Advanced non-wearing magnetic coil technology, quieter and more stable operation, longer service life.
3. Can use mobile power supply, convenient for outdoor use, not limited by experimental conditions.
4. Counterclockwise and clockwise-counterclockwise stirring functions can be switched to make the stirring more uniform and meet a variety of experimental needs.
5. The LCD screen displays the parameters in real time, the LED indicator shows the working status, and the instrument running status is more intuitive.
6. Touch buttons and sliding touch adjustment, durable, novel and stylish.
7. Dust and water proof, high protection level (IP65).

Parameters

Model	FMS-15
Max. stirring amount(H ₂ O)	1L
Motor output power	3W
Speed adjustable range	20 ~ 1500rpm
Clockwise/counterclockwise stirring intermittent time (except single direction stirring)	45s/5s
Timing Range	0 - 9hours (0 is forever)
Display method	LCD
Platform material	PVC
Platform Size	Φ110
Dimensions	W.110xD.180xH.14.5mm
weight	0.18Kgs
Voltage	DC 5V 2A
IP Rating	IP63

