
MINIP-25 Mini Magnetic Stirrer

Magnetic Series

Introduction

Mini Magnetic Stirrer is a mini-sized and portable magnetic stirrer with high standard reliability and durability. It has wide range of application.

Features

1. Distinct designed. Small and exquisite. Simple
2. operation.
3. Long using life with no carbon brush motor.
4. Speed regulation stepless, large range of speed.
5. Strong magnetism to avoid stirrer jumping.
Convenient and portable. Meet different kinds of experimental demands.

Parameters

Model	MINIP-25
Platform Size	Φ88mm
Platform Material	Plastic
Speed Range	0~2500 rpm
Max. Capacity	250ml
Stirring bar Dimension	Φ8x46mm
Motor Input/Output	1.2W/0.8W
Voltage	DC12V 1A
Power	≤12W
Dimension	W.110 x D.110 x H.70mm
Net Weight	1.2kg

